


SOLEIL status update

Martin Savko

savko@synchrotron-soleil.fr



- energy: 2.75 GeV
 - current: 500 mA
 - electron beam lifetime: ~11 hours
 - circumference: 354 m
 - emittance (horizontal, vertical): 3.7×10^{-9} , 11×10^{-12} m.rad
 - brilliance: 10^{20} ph.s⁻¹.mrad⁻¹.mm⁻² @ 0.1% bandwidth
-
- founded in 2001, in operation since 2006
 - funded jointly by CNRS (72%) and CEA (28%)
 - 350 employees



Proxima 1

Source: **U20** in vacuum undulator

Focussing: KB, **CRL**, **20x40 μm** , **project for new KB mirrors**

Tunable: Si 111 CCM, 5.5 - 15.5 keV

Flux: **2.0e12 ph/s** @ 500mA @ 12.65keV

Area Detector: **Eiger X 16M**

XRF Detector: Ketek AXAS-M2 **H150** (XIA)

OAV Camera: Prosilica GC 1350 (4.65 μm , 1360x1024)

Goniometer: **SmarGon**

Sample Changer: CATS (**48 cryo**, **16 ambient**) **Moving towards solution from SLS: Tell with a large dewar.**

MXCuBE: Qt4 v 2.3 (**CentOS 7**), HardwareRepository, Python 2.7. **mxcube**core, **mxcube**web in development (nearly ready)

Proxima 2A

Source: **U24** in vacuum undulator

Focussing: KB, **horizontal PFM**, **5x10 μm**

Tunable: Si 111 CCM, 5.5 - 18.5 keV

Flux: **1.6e12 ph/s** @ 500mA @ 12.65keV

Area Detector: **Eiger X 9M**

XRF Detector: Ketek AXAS-M2 **H80** (XIA, Xpress3)

OAV Camera: **BZoom**

Goniometer: **MD3Up with minikappa (MK3 and HC/REX)**, **Plate Screener**

Sample Changer: CATS (**144 cryo**, **48 ambient**)

MXCuBE: Qt5 (**Ubuntu 20.04**), **mxcube**core, Python 3.8, **Adating mxcube**qt to **mxcube**core with **Abstract**Diffractionmeter, GPhL workflows, Will investigate Python $\geq$ 3.11

work done over the past 6 months

- **unattended data collection**


- available on both beamlines (Murko, fast X-ray centring, volume reconstruction)
- full ISPyB integration
- speeding up the collects down to 2'15" per sample on Proxima2
 - mount: 40 s; opti: 25 s; tomo: 40s; char 60 s; main: ~30 s per sweep
- groundwork for GPhL workflows integration into unattended mode

- **mxcubeqt**

- PX2 mxcubecore with GenericDiffractometer, moving towards AbstractDiffractometer
- PX1 based on HardwareRepository

- **mxweb**

- PX2 mxcubecore@AbstractDiffractometer mxweb@
- PX1 mxcubecore@GenericDiffractometer mxweb@
 - Used internally, 2FA authentication implemented



File Queue View Graphics Help

Sample alignment

 ω : 135.00 λ : 0.00 ϕ : 0.00

Zoom: 1


Centre Save Orient Beam Anneal XC...re

Frontlight: 35

Backlight: 0

Transfer

oav

☐ Display beam size X: 873 Y: 61
☐ Graphics items

Sample: 9_11

Standard Collection

Characterisation

Helical Collection

Energy Scan

XRF Spectrum

Data location

Folder:

/nfs/data4/2025_Run5/com-proxima2a/2025-11-18/RAW_DATA

PX2-0046

File name: pos11_1_#####.raw

Prefix

Run number 1

Parameters

Count time (s): 10

Excitation energy (keV): 15.00

☒ Adjust transmission

GPhL Workflows

Advanced

Collect Now

Add to queue

[2025-11-18 23:49:51] Collection finished
[2025-11-18 23:49:51] Collection: Updating data
collection in LIMS
[2025-11-18 23:49:51] Collection: Updating data
collection in LIMS
[2025-11-18 23:49:51] Collection finished
[2025-11-18 23:49:51] Characterisation: Please wait ...
[2025-11-18 23:49:51] Characterising, please wait ...
[2025-11-18 23:49:51] Characterisation completed.

ISPyB proposal

Logout

Sample tree

Mode: Cats

Show robot menu

Sample:

ISPyB

Centring: Auto + n-click

n-clicks: 3 step 120.0

Filter:

No filter

9_5
9_6
9_7
9_8
9_9
9_10
9_11
Characterisation - 1
refnos11_1 (Point - not defined)

☐ Queue history

Collect Queue Pause

Frontend shutter

closed

Open

Close

Safety shutter

disabled

Open

Close

Machine current

501.0 mA

Machine state

Wed Nov 19 08:06

Shift Lignes

filling: 4/4

Derniere perte : Default

Profibus Alim QPole

Beam usable

Hutch temperature

21.8 C

Flux

5.62e+10 ph/s

Beam size

0.010x0.005 mm

Cryostream

In place

sample temperature: 100.0 K

Sample changer

Dewar level in range

refill On

Ramdisk

Total: 159.4TB

Free: 3.7TB (2%)

State: Ready Diffractometer: Ready Sample changer: Last collect: Characterisation : Successful (2025-11-18 23:49:51)

File system EDNA ISPyB

FileQueueGraphicsViewHelp

CollectLogTest

ISPyB proposal

idtest0 - operator on IDTESTeh1

Sample tree

Mode: Sample changer

Show SC-details

Sample:

ISPyB

Centring: Manual 3-click

Filter:

No filter

Puck 1

1:1

1:2

1:3

1:4

1:5

1:6

1:7

1:8

1:9

1:10

Puck 2

2:1

2:2

2:3

2:4

2:5

2:6

2:7

2:8

2:9

2:10

Puck 3

3:1

3:2

3:3

3:4

3:5

Queue history

Collect Queue

Sample centring

Omega: 281.53

Phi: 6.69

Kappa: 1.16

Zoom 1

Sample video

Centre

Save

Line

Grid

Focus

Snapshot

Refresh

Align

Select all

Clear all

Auto

Display beam size X: 55 Y: 80

Graphics items

Centring Point 1 (kappa: 1.16 phi: 6.69) created

Beam size

Horizontal: 50 µm

Vertical: 50 µm

Slits

Horizontal: 50 µm

Vertical: 50 µm

Aperture

Phase

CENTRE

Standard Collection

Sample: manually-mounted

Acquisition

Oscillation start: 281.53

Osc. range per frame: 0.1

Number of images: 1

Total osc. range: 0.10

First image: 1

Full osc. range

Exposure time (s): 0.04

Total exp. time (s): 0.04

Kappa: 1.1618

Phi: 6.6856

Energy (keV): 12.4

MAD

Resolution (Å): 4.402

Detector mode: 0

Transmission (%): 99.9

Shutterless

Data location

Folder:

Characterisation

Helical Collection

Energy Scan

XRF Spectrum

GEL Workflows

Advanced

SSX

Comments

Collect Now

Add to queue

Energy

Current: 12.4000 keV

1.000 Å

Set to:

Center beam after energy change

Transmission

Current: 99.90 %

Set to:

Resolution

Current: 4.402 Å

500.00 mm

Set to:

Door interlock

unknown

Unlock

Shutter

Closed

Open

Close

Online processing results

0.0

0.1

0.2

0.3

0.4

0.5

0.6

0.7

0.8

0.9

Result: Resolution

Auto scale

Display image in ADxV

Image: #, value #

Threshold:

Create centring points

[2026-05-20 19:34:00] ISPyB proposal: idtest0 - operator on IDTESTeh1

[2026-05-20 19:34:00] Data collect is disabled

[2026-05-20 19:34:00] Safety shutter is closed (Open the safety shutter to enable collections)

[2026-05-20 19:34:07] Centring click at x:349, y:221

[2026-05-20 19:34:08] Centring click at x:314, y:229

[2026-05-20 19:34:09] Centring click at x:231, y:231

[2026-05-20 19:34:12] Centring successful

[2026-05-20 19:34:21] Centring successful

[2026-05-20 19:34:40] Select two centred position (CTRL click) to continue

[2026-05-20 19:34:44] Automatic loop centring failed

[2026-05-20 19:34:46] Centring click at x:396, y:232

[2026-05-20 19:34:47] Centring click at x:364, y:233

[2026-05-20 19:34:48] Centring click at x:294, y:243

[2026-05-20 19:34:50] Centring successful

com-proxima2a@unknownState: ReadyDiffractometer: ReadySample changer: - Last collect: -


File system EDNA ISPyB

- Excellent performance of the new goniometer
 - FATs Dec 2024; SATs Jan 2025
 - SOC ~ 150nm
 - MK3 Translation calibration < 4 μm per axis
 - MK3 no longer necessary for sample exchange
 - was 180 deg, now 30 deg operating at this angle by default
 - Cryostream repositioned
 - Nearly instantaneous zoom changes.
 - Slipring issue December 2025


MK3 translation calibration

- Expected errors:


- median errors: AlignmentY 1.7 μm , CentringX 3.0 μm , CentringY 3.7 μm
- 3d error: 6.9 μm





calibration: , optimization method: nelder-mead



Commissioning new goniometer!










8_6_3, sampling 17.1 Hz



- 

Acknowledgements

- **SOLEIL team:** Bill Shepard, Serena Sirigu, Rémi Soen, Damien Jeangerard, Eric Larquet, Pierre Legrand, Tatiana Isabet, Robin Lener, Andrew Thompson, Dan Costin, Arthur Felisaz-Calvino
- **GPhL team:** Rasmus Fogh, Peter Keller, Clemens Vonnrhein, Claus Flensburg and Gérard Bricogne
- **SLS CF** May Sharpe and Dennis Stagman
- **EMBL HH team:** Marina Novikova and Gleb Bourenkov
- **Murko collaboration:** Andrey Gruzinov and Thomas White (DESY), Didier Nurizzo (ESRF), Arthur Felisaz-Calvino (SOLEIL); David Aragao and Ralf Fleig (DLS), Annie Heroux & Nicola Demitru (Elettra), Jie Nan and Isak Lindé (Max IV), Roeland Boer (Alba), Tom Crosskey (Bessy), Kate Smith and Ezequiel Panepucci (formerly SLS), Scott Classen (ALS)