


SOLEIL Status Report


Author: Mohsen Dahesh
presented by: Idrissou Chado

Py-ISPyB development status @ SOLEIL




Task	Status	Comments
SDLC automation	80%	Branching strategies, merge request rules, code review, <u>CI/CD pipelines</u> , quality checks, and approval flows
Authentication / Authorization by Keycloak	100%	
Resolve bugs	50%	Notice bug and treating them
Shippings import csv/excel	100%	Import csv/excel files of shipment data and tested by key-user.
Shippings Container	95%	Create or edit container's samples

SDLC Automation



Shippings Import CSV or Excel files



The screenshot displays the SOLEIL Shippings management interface. The top navigation bar includes links for Sessions, Calendar, Contacts, Shipments, Proteins, and Samples. The main content area shows a shipment detail view for 'test2'. A modal dialog titled 'Upload your CSV / Excel' is open, featuring a file selection button 'Choisir un fichier', a status 'Aucun fichier choisi', and 'CANCEL' and 'UPLOAD' buttons. A red arrow points to the 'Import from CSV / Excel' button in the bottom section of the shipment detail view.

mx20240026 Sessions Calendar Contacts Shipments Proteins Samples

Home » Proposals » mx20240026 » Shipments

Shipments Shipping addresses

+ Create new shipment

search in the last 4 shipments

test2 opened
2025-03-28
0 parcels (30 samples)

Vernalis-Oct2024 processing
2024-09-25
1 parcels (80 samples)

Vernalis-June2024 at_SOLEIL
2024-06-04
1 parcels (112 samples)

Vernalis-April2024 at_SOLEIL
2024-03-27
1 parcels (112 samples)

test2
Beamline: Sender address: labContact1
Date: Return address: ROBLIN-Synchrotron SOL
Status: opened
Comments: none

One of your labels is not printed

Send shipment to the facility Edit Delete

Content (0 parcels - 30 samples) Content Statistics

+ Create parcel Import from CSV / Excel Export to pdf

Upload your CSV / Excel

Choisir un fichier Aucun fichier choisi

CANCEL UPLOAD

After importing a CSV file

Edit

 Send shipment to the facility

 Edit

 Delete

Content (0 parcels - 30 samples)

ContentStatistics

+ Create parcel

 Import from CSV / Excel


 Export to pdf

Name: Dewar1

Status:


Location:

Storage:




code: MX-001

type: Unipuck




code: MX-002

type: Unipuck




code: MX-003

type: Unipuck



code: MX-004

type: Unipuck




code: MX-005

type: Unipuck



code: MX-006

type: Unipuck



code: MX-007

type: Unipuck

Add container... | v

 Edit

 [Print Label](#)

 [Print Label - modal](#)

 Remove

Edit container's sample

Edit container

Editing container from:
test2 (Dewar1)

Code: *

Type: Unipuck

Beamline :

#Sample Changer:

Status:

⚠ Sample name field is mandatory and no special characters are allowed

⚠ Protein + sample name should be unique for the whole proposal

#	Protein Acronym	Sample Name	Pin BarCode	Crystal Form	Exp. Type	Aimed resolution	Required resolution	Beam Diameter	Number of positions	Aimed multiplicity	Aimed Completeness	Forced Space G.	Radiation Sensitivity	Smiles	Tot Rot. Angle	Min Osc. Angle	Observed resolution	Comments	Ligands	Energy	Transmission range	Exposure Time	Oscillation Range
1	bo	SampleA01		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
2	bo	SampleA02		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
3	bo	SampleA03		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
4	bo	SampleA04		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
5	bo	SampleA05		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
6	bo	SampleA06		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
7	bo	SampleA07		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
8	bo	SampleA08		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
9	bo	SampleA09		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
10	bo	SampleA050		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
11	bo	SampleA011		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
12	bo	SampleA012		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
13	bo	SampleA013		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
14	bo	SampleA014		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
15	bo	SampleA015		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1
16	bo	SampleA016		undefined - (null, null, null - null, null, null)	Default	1.5		50							360					17	100	0.002	0.1

Save


Create samples for container

Name: Dewar178

Status:

Location:

Storage: Dry Shipper



code:

type: Unipuck


Create samples for a container

Name: Dewar178

Status:

Location:

Storage: Dry Shipper



code: abcd

type: Unipuck

- **SUN set synchronization**
 - Note: SUN set is the Digital User Office tool at SOLEIL
- **Prepare experiment**
- **CryoEM**
- **Improve API performance**
 - n+1 query problem, caching, heavy computation APIs, etc
- **Update third party libraries**
- **Migration of BioSAXS features**

Thank you

And Special Thanks :

- *To A. NARDELLA, F. CORONIS, A. GONZALEZ and J. NAN **from MAXIV***
- *To OUNSY, T. ISABET, P. LEGRAND, A. THUREAU, J. PEREZ, M. SAVKO and W. SHEPARD **from SOLEIL***

